

Paed Emergency Trolley Checklist

										Month _____	Year _____	Location _____												
(Defib type: Zoll AED 3)	Time <i>(MUST be day shift)</i>	Whole Trolley Clean	Portable suction checked/ready to use	Defibrillator inspected, Green tick present on indicator	1x sealed set pads (connected to device)	Oxygen >3/4 full	Valve left ON/Flowmeter to ZERO	Top of trolley contents all present as detailed on the Paediatric Emergency Trolley Equipment List	Paediatric Basic Airway Tray	Advanced Airway Accessories	Tracheal Tube Uncuffed	Tracheal Tube cuffed	Paediatric Difficult Airway	Paediatric LV Access Tray	NG & Suction Tray	Blood bottles	Fluids/Drugs as listed	Emergency Drugs Tray	Misc Drawer contents present	EZIO present and sealed (if applicable)	Missing/Expired Items Replaced	Initial	Print Name / Comments <small>NOTE: annotate here if any equipment is replaced following routine check or emergency and re-sign to confirm trolley is ready for use</small>	
1																								
2																								
3																								
4																								
5																								
6																								
7																								
8																								
9																								
10																								
11																								
12																								
13																								
14																								
15																								
16																								
17																								
18																								
19																								
20																								
21																								
22																								
23																								
24																								
25																								
26																								
27																								
28																								
29																								
30																								
31																								

CHECK MUST BE PERFORMED AFTER 0730hrs DAILY

(Additional forms available through resuscitation service website, Look under 'R' on the intranet directory and follow the links to 'equipment').
Please print a new form each month (do not stockpile - forms are updated regularly)

Resuscitation Service - updated Dec 2020



← Scan on Phone for Defibrillator User Guide