

Paediatric Emergency Trolley Checklist

Month _____ Year _____ Location Theatres only

(Defib type: m)	Time (MUST be day shift)	Clean, tidy and dust free	Defib checked and plugged in*	Portable Suction -set up and ready*	Oxygen >3/4 full Valve left ON/Flowmeter to ZERO	Top of trolley contents all present as detailed on the Paediatric Emergency Trolley Equipment List	Paediatric Basic Airway tray	Advanced Airway Accessories	Tracheal Tube Uncuffed	Tracheal Tube cuffed	Paediatric Difficult Airway	Paediatric I, V Access Tray	NG & Suction Tray	Blood Bottles	Fluids/Drugs & IO	Emergency Drugs Tray	Misc Drawer Contents present	Initial	Print Name / Comments <i>NOTE: annotate here if any equipment is replaced following routine check or emergency and re-sign to confirm trolley is ready for use</i>
1																			
2																			
3																			
4																			
5																			
6																			
7																			
8																			
9																			
10																			
11																			
12																			
13																			
14																			
15																			
16																			
17																			
18																			
19																			
20																			
21																			
22																			
23																			
24																			
25																			
26																			
27																			
28																			
29																			
30																			
31																			

* may be on adult trolley, sign to confirm these checks have been undertaken
 (Additional forms available through resuscitation service website, Look under 'R' on the intranet directory and follow the links to 'equipment').