

Adult Emergency Trolley Checklist (With Paed Basic Airway Tray & BVM)

Month _____ Year _____ Location _____

(Defib type: Zoll R Series)	Time (MUST be day shift)	Whole Trolley Clean	Portable suction checked/ready to use	Defibrillator checked/plugged in	1x sealed set adult pads (connected to device)	1 x spare set adult pads	2 x sets paed pads	ECG Electrodes	Spare pack ECG paper	Oxygen -3/4 full	Valve left ON/flowmeter to ZERO	BVM	Bougies/Introducer (loose item)	Spare Yankour Sucker	Suction tubing (Spare sealed pack)	Size 3 mask	Size 4 mask	Non rebreath mask	Suction Catheters	Sharps Bin	Gloves	Basic Airway Tray	IV Access Tray	Advanced Airway Tray	Blue Drugs Box	Red Drugs Box (if applicable)	Paediatric Basic Airway & BVM	Missing/Expired Items Replaced	Initial	Print Name / Comments <small>NOTE: annotate here if any equipment is replaced following routine check or emergency and re-sign to confirm trolley is ready for use</small>	
1																															
2																															
3																															
4																															
5																															
6																															
7																															
8																															
9																															
10																															
11																															
12																															
13																															
14																															
15																															
16																															
17																															
18																															
19																															
20																															
21																															
22																															
23																															
24																															
25																															
26																															
27																															
28																															
29																															
30																															
31																															

CHECK MUST BE PERFORMED AFTER 0730hrs DAILY
 (Additional forms available through resuscitation service website, Look under 'R' on the intranet directory and follow the links to 'equipment').
 Please print a new form each month (do not stockpile - forms are updated regularly)

Resuscitation Service - updated Dec 2020



← Scan on Phone for Defibrillator User Guide